

**NULINE**  
REFRIGERATION

A SIGN OF QUALITY

**NULINE**  
+ MEDICAL

# Blanket & Fluid Warming

Nuline's comfort range of Blanket & Fluid warming cabinets are designed to maintain blankets and fluid at a safe temperature to assist patient comfort in maternity, theatre, aged care etc. The fan forced system allows even temperature throughout the cabinet. The LED digital display allows you to see the temperature of your product easily. Fitted with a built in HI alarm in case of overheating. Designed and built in Australia finished with stainless interior/exterior and chrome plated shelves.



## Features & Benefits

- Digital display shows product temperature
- Chrome plated shelves
- Stainless steel interior/ exterior
- Self closing doors
- Castors

Model	Dimension	Capacity	Shelves	Amps
BW5/FW5	500W x 525D x 600H	179 Litre	3	3.9
BW10/FW10	850W x 572D x 900H	295 Litre	4	3.9
BW21/FW21	750W x 650D x 1800H	648 Litre	6	3.9

\*Nuline Refrigeration reserves the right to alter specifications without prior notice. \*Height does not include castors. \*Voltage is 240v.