



# MEDIFRIDGE

DISPLAY REFRIGERATION, FREEZERS & HEALTH CARE PRODUCT

## NLAB Medical & Laboratory



### Optional

Chart Recorder  
Data Logger  
Auto Dialer  
Intrinsically Safe

### Specification

| Model       | Dimensions           | Capacity | Shelves | Door         |
|-------------|----------------------|----------|---------|--------------|
| NLAB 700/1  | 870 x 700d x 2110h   | 5        | 700     | Glass/ Solid |
| NLAB 1000/2 | 1200w x 700d x 2110h | 5        | 1000    | Glass/ Solid |
| NLAB 1614/3 | 1800w x 700d x 2110h | 15       | 1614    | Glass/ Solid |

### Features

Audible Alarm  
Lockable Handle  
Automatic Water Disposal  
Voltage Free Alarm Connection  
Digital Temperature Indicator  
LED Illumination