



MEDIFRIDGE

DISPLAY REFRIGERATION, FREEZERS & HEALTH CARE PRODUCT

NLMB Laboratory and Medical



Optional

Chart recorder
Data logger
Auto dialler
Voltage free alarm connection
Intrinsically safe

Features

Audible alarm
Lockable handle
Automatic Water disposal
Double glazed doors
Interior light
Full Length door handles
Digital Temperature indicator
White colourbond interior/exterior
with stainless steel floor

Specification

Model	Dimensions	Capacity	Shelves	Door
NLMB/1	750 x 700 x 1980	488	5	1 Door Glass/Solid
NLMB/2	1200 x 700 x 1980	820	5	2 Door Glass/Solid
NLMB/3	1660 x 700 x 1980	1171	15	3 Door Glass/Solid